



BIOMOLECULE CHARACTERIZATION & PURIFICATION — AT — EXTREMELY AFFORDABLE — PRICES —

AGGREGATES • ADC • mAb • BIOSIMILARS

Архангельск (8182)63-90-72
Астана (7172)727-132
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Тенза (8412)22-31-16

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13

Сургут (3462)77-98-35
Тверь (4822)63-31-35
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Жабаровск (4212)92-98-04
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Ярославль (4852)69-52-93

The background of the slide is an underwater scene. At the top, there are bright, shimmering light rays and bubbles, suggesting sunlight filtering through the water. The water transitions from a bright blue at the top to a darker, deep blue towards the bottom. In the lower half, there are faint, ethereal patterns that look like light trails or perhaps the movement of small organisms or particles in the water.

Welcome to Affordable, Higher Resolution Size Exclusion Chromatography

With Yarra™ SEC Columns you can:

- **Save Money** with extremely affordable prices
- **Achieve Better Results** through larger exclusion ranges and higher efficiencies
- **Get Enhanced Recovery** using more inert Yarra particles and Bio-Inert hardware
- **Gain Time** with faster more productive HPLC/UHPLC runs on the new 1.8 μm SEC-X150 and X300
- **Feel at Ease** knowing you have an unmatched product guarantee

Better Protein/Peptide SEC Starts With Yarra

Affordable Pricing and Unmatched Product Guarantee 4
Dependable Results 5

NEW 1.8 μm with Bio-Inert Hardware	Exceptional Separation Power (Insulin / mAbs)	6-7
	Expanded Separation Window	8
	Better Recovery with Greater Inertness	9
	HPLC/UHPLC System Versatility	10
	Incredible Performance Upgrade	11
	Excellent Lifetime and Stability	12
3 μm	Simple Upgrade for 5 μm SEC Methods	11
	Excellent Lifetime and Stability	12
5 μm	High Performance Preparative SEC	13

Select the Right Column for Your Application 14
Ordering Information 15

Affordable Pricing and Unmatched Product Guarantee

Yarra SEC columns are highly affordable and are accompanied by a product guarantee that ensures performance satisfaction, making it extremely easy to upgrade current SEC methods!

guarantee

If Yarra™ analytical columns do not provide you with at least an equivalent separation as any other GFC column of similar porosity, type, and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

Yarra vs. Waters® BEH SEC



SEC-X150		SEC-X300		VS.	Waters* 1.7 µm 150 x 4.6 mm	
					BEH125 SEC	BEH200 SEC
1.8		1.8		Particle Size (µm)	1.7	1.7
150		300		Pore Size (Å)	125	200
1K-450K		10K-700K		MW Range in Native Conditions (Da)	1K-80K	10K-450K
>30,000		>30,000		Efficiency (plates/column)	>30,000	>30,000

*Waters specifications taken from Waters website.

Yarra vs. Tosoh Bioscience® TSKgel®



Yarra 300 x 7.8 mm			VS.	TSKgel® 300 x 7.8 mm		
SEC-2000	SEC-3000	SEC-4000		G2000SWxl	G3000SWxl	G4000SWxl
3	3	3	Particle Size (µm)	5	5	8
145	290	500	Pore Size (Å)	125	250	450
1K-300K	5K-700K	15K-1,500K	MW Range in Native Conditions (Da)	5K-150K	10K-500K	20K-7,000K
48,000	48,000	38,000	Efficiency (minimum theoretical plates)	20,000	20,000	16,000

*All TSKgel specifications were taken from Tosoh Bioscience 2004-5 Laboratory Products Catalog

Comparative separations may not be representative of all applications.


Expert Manufacturing for Dependable Performance and Reliability



Reproducible GFC columns require extreme detail in every aspect of the manufacturing and packing process. First, Yarra silica particles are synthesized using narrow tolerances for pore and particle size. Next, bonding of a proprietary hydrophilic ligand is tightly controlled and packing of each column requires validated recipes tested to high specifications to ensure excellent reproducibility.

Batch-to-Batch Reproducibility


Yarra 1.8µm SEC-X150



Column: Yarra 1.8µm SEC-X150
Dimensions: 150 x 4.6 mm
Part No.: 00F-4631-E0
Mobile Phase: 100 mM Sodium Phosphate in Water
pH 6.8 + 0.025% NaN₃
Flow Rate: 0.3 mL/min
Detection: UV @ 280 nm
Sample: 1. Thyroglobulin (669 kDa)
2. IgA (300 kDa)
3. IgG (150 kDa)
4. Ovalbumin (44 kDa)
5. Myoglobin (17 kDa)
6. Uridine

Batch-to-Batch Reproducibility

Yarra 3µm SEC-3000





Column: Yarra 3µm SEC-3000
Dimensions: 300 x 7.8 mm
Part No.: 00H-4513-K0
Mobile Phase: 100 mM Sodium Phosphate in Water
pH 6.8 + 0.025% NaN₃
Flow Rate: 1 mL/min
Detection: UV @ 280 nm
Sample: 1. Thyroglobulin (669 kDa)
2. IgA (300 kDa)
3. IgG (150 kDa)
4. Ovalbumin (44 kDa)
5. Myoglobin (17 kDa)
6. Uridine

Exceptional Separation Power - Insulin

The high resolution provided by Yarra 1.8 μm SEC-X150 columns make them an excellent solution for the characterization of biotherapeutic insulin and any resulting aggregates. With two column lengths, you have the ability to further strengthen resolution with a 300 mm or speed up methods with a 150 mm.

Degraded Human Insulin




Conditions for all columns:

Column: Yarra 1.8 μm SEC-X150
BEH125 SEC, 1.7 μm
Dimensions: 150 x 4.6 mm

Mobile Phase: 1 g/L L-Arginine in Water/Acetic Acid/Acetonitrile (65:15:20)

Flow Rate: 0.2 mL/min
Detection: UV @ 280 nm
Sample: Human Insulin

Expanded Resolution of Recombinant Human Insulin



Column: Yarra 1.8 μm SEC-X150
Dimensions: 300 x 4.6 mm
Part No.: 00H-4631-E0
Mobile Phase: Arginine (1 g/L) in Water/Acetic Acid/Acetonitrile (65:15:20)
Flow Rate: 0.4 mL/min

Temperature: Ambient
Detection: UV @ 276 nm
Sample: 1. Insulin Dimer
2. Insulin Monomer


Comparative separations may not be representative of all applications.

Exceptional Separation Power— Monoclonal Antibodies (mAb)



The large separation range of the Yarra SEC-X300 makes it ideal for high molecular weight aggregation studies of monoclonal antibodies. While the smaller MW focus of the SEC-X150 allows for excellent identification and separation of a mAb and its variable region fragments.


Mouse mAb (monoclonal antibody)



Column: Yarra 1.8 μ m SEC-X300
Dimensions: 150 x 4.6 mm
Part No.: 00F-4743-E0
Mobile Phase: 0.1 M Sodium phosphate in Water (pH 6.8)
Flow Rate: 0.35 mL/min

Temperature: Ambient
Detection: UV @ 280 nm
Sample: Mouse Monoclonal antibody

Monoclonal Antibody and Fragments




Column: Yarra 1.8 μ m SEC-X150
Dimensions: 300 x 4.6 mm
Part No.: 00H-4631-E0
Mobile Phase: 0.2 M Sodium phosphate in Water (pH 6.8)
Flow Rate: 0.2 mL/min

Temperature: Ambient
Detection: UV @ 280 nm
Sample: 1. mAb(monomer)
2. Fragment 1
3. Fragment 2

Expanded Separation Window


An expanded exclusion range was expertly developed and fine-tuned for both Yarra™ 1.8 µm products. Together they offer a wider applicability for popular biomolecule separations. The SEC-X150 focuses on smaller MW separations in the range of 1 kDa-450 kDa while the SEC-X300 targets those high MW compound mixes inside of 10 kDa-700 kDa.



Low MW Separation Mix

Column: Yarra 1.8 µm SEC-X150
Dimensions: 150 x 4.6 mm
Part No.: 00F-4631-E0
Mobile Phase: 100 mM Sodium Phosphate in Water pH 6.8 + 0.03% Na₂S
Flow Rate: 0.35 mL/min
Temperature: Ambient
System: Agilent® 1260 UHPLC
Detection: UV @ 280 nm
Sample: 1. BSA Dimer (132 kDa)
2. BSA (66 kDa)
3. Ovalbumin (44 kDa)
4. Myoglobin Dimer (34 kDa)
5. Myoglobin (17 kDa)
6. Ribonuclease A (13.7 kDa)
7. Insulin (5800 Da)
8. Uridine (244 Da)

mAb with Dimer



Column: Yarra 1.8 µm SEC-X300
Dimensions: 150 x 4.6 mm
Part No.: 00F-4743-E0
Mobile Phase: 0.1 M Potassium phosphate in Water (pH 6.8)
Flow Rate: 0.35 mL/min


Temperature: Ambient
Detection: UV @ 280 nm
Sample: 1. Dimer
2. Monomer

Better Recovery with Greater Inertness




Phenomenex's proprietary surface chemistry provides an inertness that is hard to match by other GFC columns. The result is minimal adsorption of peptides, proteins, and other protein aggregates leading to more accurate quantitation.

Yarra® 1.8µm SEC-X150



Waters® BEH200 SEC 1.7µm



Conditions for all columns:

- Columns:** Yarra 1.8µm SEC-X150
BEH200 SEC 1.7µm
- Dimensions:** 150 x 4.6 mm
- Mobile Phase:** 50 mM Sodium Phosphate in Water pH 5.0/
0.1 M Sodium Sulfate
- Flow Rate:** 0.3 mL/min
- Temperature:** Ambient
- Detection:** UV @ 280 nm
- Sample:**
 1. IgG Dimer
 2. IgG Monomer
 3. Albumin



New Bio-Inert Hardware

In addition to a highly inert particle chemistry, Yarra 1.8µm SEC-X150 columns utilize a Bio-Inert hardware to further ensure accurate recoveries.

Fast GFC/SEC Methods on Your HPLC or UHPLC

Utilize the high performance of the new Yarra™ 1.8µm on the system(s) that you know and are comfortable with! The amazing separation power of the Yarra 1.8µm can easily be utilized on any HPLC and UHPLC system.



Waters® ACQUITY® UPLC®



Shimadzu® Nexera®



Agilent® 1200



Questions About LC System Compatibility?


Contact your local Phenomenex
representative for guidance!

Performance Upgrade




Upgrading to better performance from existing SEC columns is easy with the Yarra column line. Yarra 3 μm particle columns are a quick replacement for existing 5 μm columns, giving increased performance and better peak shape without method adjustment. However, large gains in analyses times, resolution and efficiency can be realized with the Yarra 1.8 μm SEC-X150 on the desired HPLC/UHPLC system.


Conventional 5 μm 2000 300 x 7.8 mm



Yarra 3 μm SEC-2000 300 x 7.8 mm



Yarra 1.8 μm SEC-X150 150 x 4.6 mm



Conditions for all columns:

Columns: Conventional 5 μm 2000
Yarra 3 μm SEC-2000
Yarra 1.8 μm SEC-X150

Dimensions: As noted in chromatograms

Mobile Phase: 100 mM Sodium Phosphate in Water pH 6.8 + 0.025% NaN_3

Flow Rate: 1 mL/min

0.3 mL/min (Yarra 1.8 μm SEC-X150 150 x 4.6 mm)

Temperature: Ambient

Detection: UV @ 280 nm

Sample: 1. Thyroglobulin (669 kDa)

2. IgA (300 kDa)

3. IgG (150 kDa)

4. Ovalbumin (44 kDa)

5. Myoglobin (17 kDa)


6. Uridine

Comparative separations may not be representative of all applications.

Excellent Stability and Lifetime

The great care exhibited in both engineering and packing of Yarra™ columns can produce stable performance over longer lifetimes than other size exclusion columns on the market. By gaining extended lifetime with Yarra columns, you save significant time and money for your lab.


Increased Lifetime Under Extreme Conditions (10% IPA and higher flow rate)



Column: Yarra 1.8 µm SEC-X150
Dimension: 150 x 4.6 mm
Part No.: 00F-4631-E0
Mobile Phase: 100 mM Sodium Phosphate in Water pH 6.8 + 0.025% Na₃ + 10% IPA
Flow Rate: 0.4 mL/min
Temperature: Ambient
Detection: UV @ 280 nm
Sample: 1. Thyroglobulin (669 kDa)
 2. IgA (300 kDa)
 3. IgG (150 kDa)
 4. Ovalbumin (44 kDa)
 5. Myoglobin (17 kDa)
 6. Uridine

Minimal change in retention and performance

Exceptional Stability Over 900 Injections





Column: Yarra 3 µm SEC-2000
Dimensions: 300 x 7.8 mm
Part No.: 00H-4512-K0
Guard Cartridge: AJ0-4487
Guard Holder: KJO-4282, SecurityGuard Guard Cartridge Kit
Mobile Phase: 100 mM Sodium Phosphate in Water pH 6.8 + 0.025% Na₃
Flow Rate: 1 mL/min
Detection: UV @ 280 nm
Sample: 1. Thyroglobulin (669 kDa)
 2. IgA (300 kDa)
 3. IgG (150 kDa)
 4. Ovalbumin (44 kDa)
 5. Myoglobin (17 kDa)
 6. Uridine

High Performance Preparative SEC



Enjoy the same selectivity and ultra-high efficiency of Yarra 3 μm for your preparative gel filtration applications. Yarra SEC PREP features a 5 μm particle size version of each Yarra 3 μm chemistry in a 21.2 mm ID column for preparative purification, desalting, and characterization of biomolecules. Yarra 5 μm PREP is available at an affordable price while maintaining the high performance given with the analytical columns.



Conditions for both applications:



Columns: Yarra 5 μm SEC-2000 PREP
Dimensions: 300 x 21.2 mm
Part No.: 00H-4619-PO
Mobile Phase: 100 mM Sodium Phosphate in Water (pH 6.8)

Flow Rate: 5 mL/min
Temperature: Ambient
Detection: UV @ 280 nm (ambient)
Sample: 1. Bovine Serum Albumin (66 kDa)
2. Myoglobin (17 kDa)

Select the Right Yarra Column for Your Application

With five Yarra™ phases covering a large exclusion window, the Yarra product line offers the perfect solution for HPLC/UHPLC size exclusion chromatography of biomolecules.

Molecular Weight (MW) Separation Ranges




Looking for a Yarra 5 µm Analytical Column?

We recommend using a 3 µm column to improve the resolution/sensitivity of your 5 µm method.

If method optimization is not an option, our BioSep™ 5 µm GFC/SEC columns are a guaranteed alternative to your current 5 µm column.

Visit www.phenomenex.com/biosep or contact your local representative for more details.



Ordering Information



Yarra Column Specifications

	SEC-X150	SEC-X300	SEC-2000	SEC-3000	SEC-4000	SEC-2000 PREP	SEC-3000 PREP	SEC-4000 PREP
Resin Type	Silica	Silica	Silica	Silica	Silica	Silica	Silica	Silica
Particle Size (µm)	1.8	1.8	3	3	3	5	5	5
Pore Size (Å)	150	300	145	290	500	145	290	500
MW Range (Da)	1K-450K	10K-700K	1K-300K	5K-700K	15K-1500K	1K-300K	5K-700K	15K-1500K
pH Range	1.5 - 8.5	1.5 - 8.5	2.5 - 7.5	2.5 - 7.5	2.5 - 7.5	2.5 - 7.5	2.5 - 7.5	2.5 - 7.5
Typical Backpressure (psi)	3,250	3,250	1,300	1,300	1,000	800	800	800
Max Backpressure (psi)	7,000	7,000	3,000	3,000	1,700	1,500	1,500	1,500
Typical Efficiency	> 30,000 (150 x 4.6 mm)	> 30,000 (150 x 4.6 mm)	48,000 (300 x 7.8 mm)	48,000 (300 x 7.8 mm)	38,000 (300 x 7.8 mm)	30,000 (300 x 21.2 mm)	30,000 (300 x 21.2 mm)	30,000 (300 x 21.2 mm)
Max. Flow Rate (mL/min)	0.4	0.4	1.5	1.5	1.2	10	10	10

Yarra 1.8 µm SEC Bio-Inert Columns (mm)

	Analytical	Analytical
Phases	150 x 4.6 mm	300 x 4.6 mm

Yarra 1.8 µm SEC-X150	00F-4631-E0	00H-4631-E0
Yarra 1.8 µm SEC-X300	00F-4743-E0	00H-4743-E0



guarantee

If Yarra® analytical columns do not provide you with at least an equivalent separation as any other GFC column of similar porosity, type, and dimensions, return the column with comparative data within 45 days for a FULL REFUND.

Yarra 1.8 µm SEC Stainless Steel Columns (mm)

	Analytical	Analytical	SecurityGuard ULTRA Cartridges***
Phases	150 x 4.6	300 x 4.6	3/pk
Yarra 1.8 µm SEC-X150	00F-4631-E0-SS	00H-4631-E0-SS	AJ0-9512
Yarra 1.8 µm SEC-X300	00F-4743-E0-SS	00H-4743-E0-SS	AJ0-9513

For Stainless Steel Only

***SecurityGuard ULTRA cartridges require holder, Part No.: AJ0-9000

Yarra 3 µm SEC Columns (mm)

	Narrow Bore	Analytical	Analytical	SecurityGuard™ Cartridges (mm)
Phases	300 x 4.6	150 x 7.8	300 x 7.8	4 x 3.0* /10 pk
Yarra 3 µm SEC-2000	00H-4512-E0	00F-4512-K0	00H-4512-K0	AJ0-4487
Yarra 3 µm SEC-3000	00H-4513-E0	00F-4513-K0	00H-4513-K0	AJ0-4488
Yarra 3 µm SEC-4000	00H-4514-E0	—	00H-4514-K0	AJ0-4489

for ID 4.6–7.8 mm

*SecurityGuard Analytical cartridges require holder, Part No.: KJ0-4282

Yarra 5 µm PREP SEC Columns (mm)

	Preparative	SecurityGuard™ Cartridges (mm)
Phases	300 x 21.2	15 x 21.2** /ea
Yarra 5 µm SEC-2000 PREP	00H-4619-PO	AJ0-8588
Yarra 5 µm SEC-3000 PREP	00H-4620-PO	AJ0-8589
Yarra 5 µm SEC-4000 PREP	00H-4621-PO	AJ0-8950

**PREP SecurityGuard™ Cartridges require holder, Part No.: AJ0-8223 for ID 18–29 mm



Архангельск (8182)63-90-72
Астана (7172)727-132
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13

Сургут (3462)77-98-35
Тверь (4822)63-31-35
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Ярославль (4852)69-52-93