

Архангельск (8182)63-90-72
 Астана (7172)727-132
 Астрахань (8512)99-46-04
 Барнаул (3852)73-04-60
 Белгород (4722)40-23-64
 Брянск (4832)59-03-52
 Владивосток (423)249-28-31
 Волгоград (844)278-03-48
 Вологда (8172)26-41-59
 Воронеж (473)204-51-73
 Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
 Ижевск (3412)26-03-58
 Казань (843)206-01-48
 Калининград (4012)72-03-81
 Калуга (4842)92-23-67
 Кемерово (3842)65-04-62
 Киров (8332)68-02-04
 Краснодар (861)203-40-90
 Красноярск (391)204-63-61
 Курск (4712)77-13-04
 Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
 Москва (495)268-04-70
 Мурманск (8152)59-64-93
 Набережные Челны (8552)20-53-41
 Нижний Новгород (831)429-08-12
 Новокузнецк (3843)20-46-81
 Новосибирск (383)227-86-73
 Омск (3812)21-46-40
 Орел (4862)44-53-42
 Оренбург (3532)37-68-04
 Пенза (8412)22-31-16

Пермь (342)205-81-47
 Ростов-на-Дону (863)308-18-15
 Рязань (4912)46-61-64
 Самара (846)206-03-16
 Санкт-Петербург (812)309-46-40
 Саратов (845)249-38-78
 Севастополь (8692)22-31-93
 Симферополь (3652)67-13-56
 Смоленск (4812)29-41-54
 Сочи (862)225-72-31
 Ставрополь (8652)20-65-13

Сургут (3462)77-98-35
 Тверь (4822)63-31-35
 Томск (3822)98-41-53
 Тула (4872)74-02-29
 Тюмень (3452)66-21-18
 Ульяновск (8422)24-23-59
 Уфа (347)229-48-12
 Хабаровск (4212)92-98-04
 Челябинск (351)202-03-61
 Череповец (8202)49-02-64
 Ярославль (4852)69-52-93

www.phenomenex.nt-rt.ru || pxp@nt-rt.ru

Колонки капиллярные Zebron GC



Essentials

Smart starting point for the GC method developer.

ZB-1
 ZB-1ms
 ZB-5
 ZB-5ms
 ZB-35
 ZB-50
 ZB-624
 ZB-1701
 ZB-WAX
 ZB-FFAP
 ZB-XLB



Plus

Dare to defy high boilers, contaminants, and carry-overs.

ZB-1^{PLUS}
 ZB-5^{PLUS}
 ZB-5MS^{PLUS}
 ZB-WAX^{PLUS}
 ZB-624^{PLUS}



Inferno™

Best upgrades compared to their Essentials counterparts.

ZB-1HT
 ZB-5HT
 ZB-35HT
 ZB-XLB-HT



Unlimited

Unleash the power of selectivity for targeted performance that breaks from the mold.

ZB-FAME
 ZB-SemiVolatiles
 ZB-MultiResidue--1 & -2
 ZB-CLPesticides-1 & -2
 ZB-Drug-1
 ZB-BAC-1 & -2
 ZB-1XT SimDist
 ZB-Bioethanol
 ZB-DHA-PONA
 ZB-PAH-EU
 ZB-PAH-CT
 ZB-Dioxin

Наименование	Фаза	Длина (м)	Внутренний диаметр (мм)	Толщина пленки (мкм)
7CK-G001-35	ZB-1	10	0.53	2.65
7EG-G001-02	ZB-1	15	0.25	0.1
7EG-G001-11	ZB-1	15	0.25	0.25
7EG-G001-22	ZB-1	15	0.25	1
7EK-G001-05	ZB-1	15	0.53	0.15
7EK-G001-17	ZB-1	15	0.53	0.5
7EK-G001-28	ZB-1	15	0.53	1.5
7EK-G001-36	ZB-1	15	0.53	3
7EK-G001-39	ZB-1	15	0.53	5
7EM-G001-11	ZB-1	15	0.32	0.25
7EM-G001-22	ZB-1	15	0.32	1
7HG-G001-02	ZB-1	30	0.25	0.1
7HG-G001-02-B	ZB-1	30	0.25	0.1
7HG-G001-11	ZB-1	30	0.25	0.25
7HG-G001-17	ZB-1	30	0.25	0.5
7HG-G001-22	ZB-1	30	0.25	1
7HG-G001-22-B	ZB-1	30	0.25	1
7HK-G001-17	ZB-1	30	0.53	0.5
7HK-G001-17-B	ZB-1	30	0.53	0.5
7HK-G001-28	ZB-1	30	0.53	1.5

Наименование	Фаза	Длина (м)	Внутренний диаметр (мм)	Толщина пленки (мкм)
7HK-G001-28-B	ZB-1	30	0.53	1.5
7HK-G001-36	ZB-1	30	0.53	3
7HK-G001-36-B	ZB-1	30	0.53	3
7HK-G001-39	ZB-1	30	0.53	5
7HK-G001-39-B	ZB-1	30	0.53	5
7HM-G001-11	ZB-1	30	0.32	0.25
7HM-G001-11-B	ZB-1	30	0.32	0.25
7HM-G001-17	ZB-1	30	0.32	0.5
7HM-G001-22	ZB-1	30	0.32	1
7HM-G001-22-B	ZB-1	30	0.32	1
7HM-G001-36	ZB-1	30	0.32	3
7HM-G001-39	ZB-1	30	0.32	5
7HM-G001-39-B	ZB-1	30	0.32	5
7JG-G001-17	ZB-1	50	0.25	0.5
7JK-G001-39	ZB-1	50	0.53	5
7JM-G001-17	ZB-1	50	0.32	0.5
7KG-G001-02	ZB-1	60	0.25	0.1
7KG-G001-11	ZB-1	60	0.25	0.25
7KG-G001-22	ZB-1	60	0.25	1
7KG-G001-22-B	ZB-1	60	0.25	1
7KK-G001-28	ZB-1	60	0.53	1.5
7KM-G001-11	ZB-1	60	0.32	0.25
7KM-G001-11-B	ZB-1	60	0.32	0.25
7KM-G001-22	ZB-1	60	0.32	1
7KM-G001-22-B	ZB-1	60	0.32	1
7KM-G001-36	ZB-1	60	0.32	3
7MG-G001-17	ZB-1	100	0.25	0.5
7EG-G006-11	ZB-1701	15	0.25	0.25
7EM-G006-11	ZB-1701	15	0.32	0.25
7EM-G006-11-B	ZB-1701	15	0.32	0.25
7FD-G006-08	ZB-1701	20	0.18	0.18
7HG-G006-11	ZB-1701	30	0.25	0.25
7HG-G006-22	ZB-1701	30	0.25	1
7HK-G006-22	ZB-1701	30	0.53	1
7HK-G006-22-B	ZB-1701	30	0.53	1
7HM-G006-11	ZB-1701	30	0.32	0.25
7HM-G006-22	ZB-1701	30	0.32	1
7HM-G006-22-B	ZB-1701	30	0.32	1
7KG-G006-11	ZB-1701	60	0.25	0.25
7KK-G006-22	ZB-1701	60	0.53	1
7KM-G006-11	ZB-1701	60	0.32	0.25
7HG-G012-11	ZB-1701P	30	0.25	0.25
7HM-G012-11	ZB-1701P	30	0.32	0
7AK-G014-02	ZB-1HT	5	0.53	0.1
7AK-G014-02-B	ZB-1HT	5	0.53	0.1
7CM-G014-11	ZB-1HT	10	0.32	0.25
7EG-G014-02	ZB-1HT	15	0.25	0.1
7EG-G014-02-B	ZB-1HT	15	0.25	0.1
7EG-G014-11	ZB-1HT	15	0.25	0.25
7EK-G014-05	ZB-1HT	15	0.53	0.15
7EM-G014-02	ZB-1HT	15	0.32	0.1
7EM-G014-11	ZB-1HT	15	0.32	0.25

Наименование	Фаза	Длина (м)	Внутренний диаметр (мм)	Толщина пленки (мкм)
7FD-G014-08	ZB-1HT	20	0.18	0.18
7HG-G014-02	ZB-1HT	30	0.25	0.1
7HG-G014-11	ZB-1HT	30	0.25	0.25
7HG-G014-11-B	ZB-1HT	30	0.25	0.25
7HK-G014-05	ZB-1HT	30	0.53	0.15
7HM-G014-02	ZB-1HT	30	0.32	0.1
7HM-G014-11	ZB-1HT	30	0.32	0.25
7HM-G014-11-B	ZB-1HT	30	0.32	0.25
7CB-G011-02	ZB-1MS	10	0.1	0.1
7CD-G011-08	ZB-1MS	10	0.18	0.18
7CD-G011-08-B	ZB-1MS	10	0.18	0.18
7DE-G011-14	ZB-1MS	12	0.2	0.33
7EG-G011-11	ZB-1MS	15	0.25	0.25
7EG-G011-11-B	ZB-1MS	15	0.25	0.25
7FD-G011-08	ZB-1MS	20	0.18	0.18
7FD-G011-22	ZB-1MS	20	0.18	1
7GE-G011-14	ZB-1MS	25	0.2	0.33
7GE-G011-14-B	ZB-1MS	25	0.2	0.33
7HG-G011-02	ZB-1MS	30	0.25	0.1
7HG-G011-11	ZB-1MS	30	0.25	0.25
7HG-G011-11-B	ZB-1MS	30	0.25	0.25
7HG-G011-17	ZB-1MS	30	0.25	0.5
7HG-G011-22	ZB-1MS	30	0.25	1
7HK-G011-22	ZB-1MS	30	0.53	1
7HM-G011-11	ZB-1MS	30	0.32	0.25
7HM-G011-22	ZB-1MS	30	0.32	1
7KG-G011-11	ZB-1MS	60	0.25	0.25
7KG-G011-22	ZB-1MS	60	0.25	1
7KM-G011-22	ZB-1MS	60	0.32	1
7YE-G011-03	ZB-1MS	17	0.2	0.11
CG0-7894	ZB-1MS	20	0.18	0.18
7AD-G031-08	ZB-1plus	5	0.18	0
7CB-G031-02	ZB-1plus	10	0.1	0.1
7DE-G031-14	ZB-1plus	12	0.2	0.33
7EG-G031-11	ZB-1plus	15	0.25	0.25
7EG-G031-11-B	ZB-1plus	15	0.25	0.25
7EM-G031-11	ZB-1plus	15	0.32	0.25
7FD-G031-08	ZB-1plus	20	0.18	0.18
7HG-G031-02	ZB-1plus	30	0.25	0.1
7HG-G031-11	ZB-1plus	30	0.25	0.25
7HM-G031-11	ZB-1plus	30	0.32	0.25
7KD-G031-08-C	ZB-1plus	60	0.18	0.18
7KG-G031-11	ZB-1plus	60	0.25	0.25
7KG-G031-22	ZB-1plus	60	0.25	1
7AK-G026-05	ZB-1XT SimDist	5	0.53	0.15
7AK-G026-49	ZB-1XT SimDist	5	0.53	0.88
7AK-G026-55	ZB-1XT SimDist	5	0.53	0.09
7AK-G026-55-B	ZB-1XT SimDist	5	0.53	0.09
7CK-G026-05	ZB-1XT SimDist	10	0.53	0.15
7CK-G026-35	ZB-1XT SimDist	10	0.53	2.65
7CK-G026-49	ZB-1XT SimDist	10	0.53	0.88
7EK-G026-11	ZB-1XT SimDist	15	0.53	0.25

Наименование	Фаза	Длина (м)	Внутренний диаметр (мм)	Толщина пленки (мкм)
7KM-G039-11	ZB-23	60	0.32	0.25
7CB-G003-02	ZB-35	10	0.1	0.1
7EG-G003-11	ZB-35	15	0.25	0.25
7EG-G003-17	ZB-35	15	0.25	0.5
7EK-G003-22	ZB-35	15	0.53	1
7HG-G003-11	ZB-35	30	0.25	0.25
7HG-G003-11-B	ZB-35	30	0.25	0.25
7HG-G003-17	ZB-35	30	0.25	0.5
7HK-G003-17	ZB-35	30	0.53	0.5
7HK-G003-22	ZB-35	30	0.53	1
7HM-G003-11	ZB-35	30	0.32	0.25
7HM-G003-17	ZB-35	30	0.32	0.5
7KG-G003-11	ZB-35	60	0.25	0.25
7KM-G003-11	ZB-35	60	0.32	0.25
7EG-G025-02	ZB-35HT	15	0.25	0.1
7EG-G025-11	ZB-35HT	15	0.25	0.25
7EM-G025-11	ZB-35HT	15	0.32	0.25
7FD-G025-08	ZB-35HT	20	0.18	0.18
7HG-G025-02	ZB-35HT	30	0.25	0.1
7HG-G025-11	ZB-35HT	30	0.25	0.25
7HG-G025-11-B	ZB-35HT	30	0.25	0.25
7HM-G025-11	ZB-35HT	30	0.32	0.25
CG0-9342	ZB-35HT	0.5	0.1	0.05
CG0-9346	ZB-35HT	0.5	0.1	0.05
7EG-G002-02	ZB-5	15	0.25	0.1
7EG-G002-11	ZB-5	15	0.25	0.25
7EG-G002-17	ZB-5	15	0.25	0.5
7EG-G002-22	ZB-5	15	0.25	1
7EK-G002-17	ZB-5	15	0.53	0.5
7EK-G002-28	ZB-5	15	0.53	1.5
7EK-G002-36	ZB-5	15	0.53	3
7EM-G002-11	ZB-5	15	0.32	0.25
7EM-G002-11-B	ZB-5	15	0.32	0.25
7EM-G002-22	ZB-5	15	0.32	1
7FD-G002-08	ZB-5	20	0.18	0.18
7FD-G002-08-B	ZB-5	20	0.18	0.18
7HG-G002-02	ZB-5	30	0.25	0.1
7HG-G002-11	ZB-5	30	0.25	0.25
7HG-G002-11-B	ZB-5	30	0.25	0.25
7HG-G002-17	ZB-5	30	0.25	0.5
7HG-G002-22	ZB-5	30	0.25	1
7HK-G002-17	ZB-5	30	0.53	0.5
7HK-G002-22	ZB-5	30	0.53	1
7HK-G002-28	ZB-5	30	0.53	1.5
7HK-G002-36	ZB-5	30	0.53	3
7HK-G002-39	ZB-5	30	0.53	5
7HK-G002-39-B	ZB-5	30	0.53	5
7HM-G002-02	ZB-5	30	0.32	0.1
7HM-G002-11	ZB-5	30	0.32	0.25
7HM-G002-11-B	ZB-5	30	0.32	0.25
7HM-G002-17	ZB-5	30	0.32	0.5
7HM-G002-22	ZB-5	30	0.32	1

Наименование	Фаза	Длина (м)	Внутренний диаметр (мм)	Толщина пленки (мкм)
7HM-G002-22-B	ZB-5	30	0.32	1
7HM-G002-36-C	ZB-5	30	0.32	3
7KG-G002-02	ZB-5	60	0.25	0
7KG-G002-11	ZB-5	60	0.25	0.25
7KG-G002-11-B	ZB-5	60	0.25	0.25
7KG-G002-17	ZB-5	60	0.25	0.5
7KG-G002-22	ZB-5	60	0.25	1
7KK-G002-28	ZB-5	60	0.53	1.5
7KM-G002-11	ZB-5	60	0.32	0.25
7KM-G002-22	ZB-5	60	0.32	1
7CB-G004-02	ZB-50	10	0.1	0.1
7CG-G004-05	ZB-50	10	0.25	0.15
7CK-G004-32	ZB-50	10	0.53	2
7EG-G004-05	ZB-50	15	0.25	0.15
7EG-G004-11	ZB-50	15	0.25	0.25
7EK-G004-22	ZB-50	15	0.53	1
7EM-G004-11	ZB-50	15	0.32	0.25
7FD-G004-08	ZB-50	20	0.18	0.18
7HG-G004-05	ZB-50	30	0.25	0.15
7HG-G004-11	ZB-50	30	0.25	0.25
7HG-G004-11-B	ZB-50	30	0.25	0.25
7HG-G004-17	ZB-50	30	0.25	0.5
7HG-G004-17-B	ZB-50	30	0.25	0.5
7HK-G004-22	ZB-50	30	0.53	1
7HM-G004-11	ZB-50	30	0.32	0.25
7HM-G004-17	ZB-50	30	0.32	0.5
7KG-G004-11	ZB-50	60	0.25	0.25
7KG-G004-17	ZB-50	60	0.25	0.5
7KK-G004-22-C	ZB-50	60	0.53	1
7PD-G004-08-C	ZB-50	40	0.18	0.18
7CM-G015-02	ZB-5HT	10	0.32	0.1
7CM-G015-02-GST	ZB-5HT	10	0.32	0.1
7EG-G015-02	ZB-5HT	15	0.25	0.1
7EG-G015-02-B	ZB-5HT	15	0.25	0.1
7EG-G015-11	ZB-5HT	15	0.25	0.25
7EK-G015-05	ZB-5HT	15	0.53	0.15
7EM-G015-02	ZB-5HT	15	0.32	0.1
7EM-G015-02-B	ZB-5HT	15	0.32	0.1
7EM-G015-02-GSA	ZB-5HT	15	0.32	0.1
7EM-G015-02-GST	ZB-5HT	15	0.32	0.1
7EM-G015-11	ZB-5HT	15	0.32	0.25
7EM-G015-11-B	ZB-5HT	15	0.32	0.25
7FD-G015-08	ZB-5HT	20	0.18	0.18
7HG-G015-02	ZB-5HT	30	0.25	0.1
7HG-G015-11	ZB-5HT	30	0.25	0.25
7HG-G015-17	ZB-5HT	30	0.25	0.5
7HK-G015-05	ZB-5HT	30	0.53	0.15
7HM-G015-02	ZB-5HT	30	0.32	0.1
7HM-G015-02-B	ZB-5HT	30	0.32	0.1
7HM-G015-11	ZB-5HT	30	0.32	0.25
7HM-G015-11-B	ZB-5HT	30	0.32	0.25
7KG-G015-11	ZB-5HT	60	0.25	0.25

Наименование	Фаза	Длина (м)	Внутренний диаметр (мм)	Толщина пленки (мкм)
7CB-G010-02	ZB-5MS	10	0.1	0.1
7CD-G010-08	ZB-5MS	10	0.18	0.18
7DE-G010-14	ZB-5MS	12	0.2	0.33
7EG-G010-11	ZB-5MS	15	0.25	0.25
7FD-G010-08	ZB-5MS	20	0.18	0.18
7FD-G010-51	ZB-5MS	20	0.18	0.32
7GE-G010-14	ZB-5MS	25	0.2	0.33
7GE-G010-14-B	ZB-5MS	25	0.2	0.33
7HG-G010-11	ZB-5MS	30	0.25	0.25
7HG-G010-11-B	ZB-5MS	30	0.25	0.25
7HG-G010-17	ZB-5MS	30	0.25	0.5
7HG-G010-22	ZB-5MS	30	0.25	1
7HG-G010-22-B	ZB-5MS	30	0.25	1
7HM-G010-11	ZB-5MS	30	0.32	0.25
7HM-G010-17	ZB-5MS	30	0.32	0.5
7HM-G010-22	ZB-5MS	30	0.32	1
7JE-G010-14	ZB-5MS	50	0.2	0.33
7KG-G010-02	ZB-5MS	60	0.25	0.1
7KG-G010-11	ZB-5MS	60	0.25	0.25
7KM-G010-11	ZB-5MS	60	0.32	0.25
7CD-G018-08	ZB-5MSi	10	0.18	0.18
7EG-G018-11	ZB-5MSi	15	0.25	0.25
7FD-G018-08	ZB-5MSi	20	0.18	0.18
7HG-G018-11	ZB-5MSi	30	0.25	0.25
7HG-G018-17	ZB-5MSi	30	0.25	0.5
7HG-G018-22	ZB-5MSi	30	0.25	1
7HM-G018-11	ZB-5MSi	30	0.32	0.25
7HM-G018-11-B	ZB-5MSi	30	0.32	0.25
7HM-G018-17	ZB-5MSi	30	0.32	0.5
7KG-G018-11	ZB-5MSi	60	0.25	0.25
7EG-G030-11	ZB-5MSplus	15	0.25	0.25
7EG-G030-17	ZB-5MSplus	15	0.25	0.5
7FD-G030-08	ZB-5MSplus	20	0.18	0.18
7FD-G030-53	ZB-5MSplus	20	0.18	0.36
7HG-G030-11	ZB-5MSplus	30	0.25	0.25
7HG-G030-11-B	ZB-5MSplus	30	0.25	0.25
7HG-G030-17	ZB-5MSplus	30	0.25	0.5
7HG-G030-17-B	ZB-5MSplus	30	0.25	0.5
7HG-G030-22	ZB-5MSplus	30	0.25	1
7HK-G030-22	ZB-5MSplus	30	0.53	1
7HM-G030-11	ZB-5MSplus	30	0.32	0.25
7HM-G030-17	ZB-5MSplus	30	0.32	0.5
7HM-G030-22	ZB-5MSplus	30	0.32	1
7HM-G030-28	ZB-5MSplus	30	0.32	1.5
7KG-G030-02-C	ZB-5MSplus	60	0.25	0.1
7KG-G030-11	ZB-5MSplus	60	0.25	0.25
7KG-G030-22	ZB-5MSplus	60	0.25	1
7KM-G030-22	ZB-5MSplus	60	0.32	1
7XG-G030-11	ZB-5MSplus	1.5	0.25	0
7CD-G032-08	ZB-5plus	10	0.18	0.18
7EG-G032-11	ZB-5plus	15	0.25	0.25
7FD-G032-08	ZB-5plus	20	0.18	0.18

Наименование	Фаза	Длина (м)	Внутренний диаметр (мм)	Толщина пленки (мкм)
7FD-G032-53	ZB-5plus	20	0.18	0.36
7HG-G032-11	ZB-5plus	30	0.25	0.25
7HG-G032-17	ZB-5plus	30	0.25	0.5
7HG-G032-17-B	ZB-5plus	30	0.25	0.5
7HG-G032-22	ZB-5plus	30	0.25	1
7HM-G032-11	ZB-5plus	30	0.32	0.25
7HM-G032-17	ZB-5plus	30	0.32	0.5
7HM-G032-22	ZB-5plus	30	0.32	1
7KG-G032-11	ZB-5plus	60	0.25	0.25
7FD-G005-22	ZB-624	20	0.18	1
7HG-G005-27	ZB-624	30	0.25	1.4
7HG-G005-27-B	ZB-624	30	0.25	1.4
7HK-G005-22	ZB-624	30	0.53	1
7HK-G005-36	ZB-624	30	0.53	3
7HK-G005-36-B	ZB-624	30	0.53	3
7HM-G005-31	ZB-624	30	0.32	1.8
7HM-G005-31-B	ZB-624	30	0.32	1.8
7KD-G005-22-C	ZB-624	60	0.18	1
7KG-G005-27	ZB-624	60	0.25	1.4
7KK-G005-36	ZB-624	60	0.53	3
7KM-G005-31	ZB-624	60	0.32	1.8
7KM-G005-31-B	ZB-624	60	0.32	1.8
7LK-G005-36	ZB-624	75	0.53	3
7NK-G005-36	ZB-624	105	0.53	3
7FD-G040-22	ZB-624plus	20	0.18	1
7FG-G040-27	ZB-624plus	20	0.25	1.4
7HG-G040-27	ZB-624plus	30	0.25	1.4
7HG-G040-27-B	ZB-624plus	30	0.25	1.4
7HK-G040-36	ZB-624plus	30	0.53	3
7HK-G040-36-B	ZB-624plus	30	0.53	3
7HM-G040-31	ZB-624plus	30	0.32	1.8
7HM-G040-31-B	ZB-624plus	30	0.32	1.8
7KG-G040-27	ZB-624plus	60	0.25	1.4
7KG-G040-27-B	ZB-624plus	60	0.25	1.4
7KK-G040-36	ZB-624plus	60	0.53	3
7KM-G040-31	ZB-624plus	60	0.32	1.8
7LK-G040-36	ZB-624plus	75	0.53	3
7KG-G037-10	ZB-88	60	0.25	0.2
7HK-G021-36	ZB-BAC-1	30	0.53	3
7HM-G021-31	ZB-BAC-1	30	0.32	1.8
7HK-G022-32	ZB-BAC-2	30	0.53	2
7HM-G022-25	ZB-BAC-2	30	0.32	1.2
7EG-G020-22	ZB-Bioethanol	15	0.25	1
7HG-G020-22	ZB-Bioethanol	30	0.25	1
7HG-G028-11	ZB-CLPesticides-1	30	0.25	0.25
7HK-G028-17	ZB-CLPesticides-1	30	0.53	0.5
7HM-G028-17	ZB-CLPesticides-1	30	0.32	0.5
7HM-G028-51	ZB-CLPesticides-1	30	0.32	0.32
KG0-9286	ZB-CLPesticides-1	30	0.53	0
KG0-9290	ZB-CLPesticides-1	30	0.3	0
7HG-G029-10	ZB-CLPesticides-2	30	0.25	0.2
7HK-G029-16	ZB-CLPesticides-2	30	0.53	0.42

Наименование	Фаза	Длина (м)	Внутренний диаметр (мм)	Толщина пленки (мкм)
7HM-G029-11	ZB-CLPesticides-2	30	0.32	0.25
7HM-G029-17	ZB-CLPesticides-2	30	0.32	0.5
7AG-G042-22	ZB-DHA-PONA	5	0.25	1
7JE-G042-17	ZB-DHA-PONA	50	0.2	0.5
7MG-G042-17	ZB-DHA-PONA	100	0.25	0.5
7QG-G042-22	ZB-DHA-PONA	150	0.25	1
7CD-G023-08	ZB-Drug-1	10	0.18	0.18
7EG-G023-11	ZB-Drug-1	15	0.25	0.25
7HG-G023-11	ZB-Drug-1	30	0.25	0.25
7FD-G033-05	ZB-FAME	20	0.18	0.15
7HG-G033-10	ZB-FAME	30	0.25	0.2
7HG-G033-10-B	ZB-FAME	30	0.25	0.2
7KG-G033-10	ZB-FAME	60	0.25	0.2
7KG-G033-10-B	ZB-FAME	60	0.25	0.2
7MG-G033-10	ZB-FAME	100	0.25	0.2
7EG-G009-11	ZB-FFAP	15	0.25	0.25
7EK-G009-22	ZB-FFAP	15	0.53	1
7EM-G009-11	ZB-FFAP	15	0.32	0.25
7HG-G009-11	ZB-FFAP	30	0.25	0.25
7HK-G009-22	ZB-FFAP	30	0.53	1
7HK-G009-22-B	ZB-FFAP	30	0.53	1
7HM-G009-11	ZB-FFAP	30	0.32	0.25
7HM-G009-17	ZB-FFAP	30	0.32	0.5
7HM-G009-17-B	ZB-FFAP	30	0.32	0.5
7HM-G009-22	ZB-FFAP	30	0.32	1
7JM-G009-17	ZB-FFAP	50	0.32	0.5
7KG-G009-11	ZB-FFAP	60	0.25	0.25
7FD-G016-08	ZB-MultiResidue-1	20	0.18	0.18
7HG-G016-11	ZB-MultiResidue-1	30	0.25	0.25
7HK-G016-17	ZB-MultiResidue-1	30	0.53	0.5
7HM-G016-11	ZB-MultiResidue-1	30	0.32	0.25
7HM-G016-17	ZB-MultiResidue-1	30	0.32	0.5
KG0-8237	ZB-MultiResidue-1	30	0.25	0
KG0-8238	ZB-MultiResidue-1	30	0.32	0
KG0-8239	ZB-MultiResidue-1	30	0.53	0
7FD-G017-08-C	ZB-MultiResidue-2	20	0.18	0
7HG-G017-10	ZB-MultiResidue-2	30	0.25	0.2
7HM-G017-11	ZB-MultiResidue-2	30	0.32	0.25
7CB-G043-59	ZB-PAH-EU	10	0.1	0.08
7FD-G043-47	ZB-PAH-EU	20	0.18	0.14
7HG-G043-10	ZB-PAH-EU	30	0.25	0.2
7KG-G043-10	ZB-PAH-EU	60	0.25	0.2
7EG-G027-11	ZB-SemiVolatiles	15	0.25	0.25
7EG-G027-17	ZB-SemiVolatiles	15	0.25	0.5
7FD-G027-08	ZB-SemiVolatiles	20	0.18	0.18
7FD-G027-53	ZB-SemiVolatiles	20	0.18	0.36
7HG-G027-11	ZB-SemiVolatiles	30	0.25	0.25
7HG-G027-11-B	ZB-SemiVolatiles	30	0.25	0.25
7HG-G027-17	ZB-SemiVolatiles	30	0.25	0.5
7HM-G027-11	ZB-SemiVolatiles	30	0.32	0.25
7KG-G027-11	ZB-SemiVolatiles	60	0.25	0.25
7KG-G027-11-B	ZB-SemiVolatiles	60	0.25	0.25

Наименование	Фаза	Длина (м)	Внутренний диаметр (мм)	Толщина пленки (мкм)
7CB-G007-02	ZB-WAX	10	0.1	0.1
7EG-G007-11	ZB-WAX	15	0.25	0.25
7EG-G007-11-B	ZB-WAX	15	0.25	0.25
7EK-G007-22	ZB-WAX	15	0.53	1
7EK-G007-22-B	ZB-WAX	15	0.53	1
7EM-G007-11	ZB-WAX	15	0.32	0.25
7FD-G007-08	ZB-WAX	20	0.18	0.18
7HG-G007-05	ZB-WAX	30	0.25	0.15
7HG-G007-11	ZB-WAX	30	0.25	0.25
7HG-G007-11-B	ZB-WAX	30	0.25	0.25
7HG-G007-17	ZB-WAX	30	0.25	0.5
7HG-G007-22	ZB-WAX	30	0.25	1
7HK-G007-17	ZB-WAX	30	0.53	0.5
7HK-G007-22	ZB-WAX	30	0.53	1
7HK-G007-22-B	ZB-WAX	30	0.53	1
7HK-G007-32-C	ZB-WAX	30	0.53	2
7HM-G007-05	ZB-WAX	30	0.32	0
7HM-G007-11	ZB-WAX	30	0.32	0.25
7HM-G007-17	ZB-WAX	30	0.32	0.5
7HM-G007-17-B	ZB-WAX	30	0.32	0.5
7JG-G007-11	ZB-WAX	50	0.25	0.25
7JM-G007-22-C	ZB-WAX	50	0.32	1
7KG-G007-05	ZB-WAX	60	0.25	0.15
7KG-G007-11	ZB-WAX	60	0.25	0.25
7KG-G007-17	ZB-WAX	60	0.25	0.5
7KK-G007-22	ZB-WAX	60	0.53	1
7KM-G007-11	ZB-WAX	60	0.32	0.25
7KM-G007-17	ZB-WAX	60	0.32	0.5
7KM-G007-17-B	ZB-WAX	60	0.32	0.5
7KM-G007-22-C	ZB-WAX	60	0.32	1
7CB-G013-02	ZB-WAXplus	10	0.1	0.1
7EG-G013-11	ZB-WAXplus	15	0.25	0.25
7EK-G013-22	ZB-WAXplus	15	0.53	1
7FD-G013-08	ZB-WAXplus	20	0.18	0.18
7HG-G013-11	ZB-WAXplus	30	0.25	0.25
7HG-G013-11-B	ZB-WAXplus	30	0.25	0.25
7HG-G013-17	ZB-WAXplus	30	0.25	0.5
7HK-G013-11	ZB-WAXplus	30	0.53	0.25
7HK-G013-11-B	ZB-WAXplus	30	0.25	0.25
7HK-G013-22	ZB-WAXplus	30	0.53	1
7HK-G013-22-B	ZB-WAXplus	30	0.53	1
7HM-G013-11	ZB-WAXplus	30	0.32	0.25
7HM-G013-11-B	ZB-WAXplus	30	0.32	0.25
7HM-G013-17	ZB-WAXplus	30	0.32	0.5
7HM-G013-17-B	ZB-WAXplus	30	0.32	0.5
7HM-G013-22	ZB-WAXplus	30	0.32	1
7HM-G013-22-B	ZB-WAXplus	30	0.32	1
7KG-G013-05	ZB-WAXplus	60	0.25	0.15
7KG-G013-11	ZB-WAXplus	60	0.25	0.25
7KG-G013-11-B	ZB-WAXplus	60	0.25	0.25
7KG-G013-17	ZB-WAXplus	60	0.25	0.5
7KK-G013-22	ZB-WAXplus	60	0.53	1

Наименование	Фаза	Длина (м)	Внутренний диаметр (мм)	Толщина пленки (мкм)
7KK-G013-22-B	ZB-WAXplus	60	0.53	1
7KM-G013-11	ZB-WAXplus	60	0.32	0.25
7KM-G013-17	ZB-WAXplus	60	0.32	0.5
7KM-G013-17-B	ZB-WAXplus	60	0.32	0.5
7FD-G019-08	ZB-XLB	20	0.18	0.18
7HG-G019-11	ZB-XLB	30	0.25	0.25
7HG-G019-17	ZB-XLB	30	0.25	0.5
7HK-G019-28	ZB-XLB	30	0.53	1.5
7HM-G019-11	ZB-XLB	30	0.32	0.25
7HM-G019-17	ZB-XLB	30	0.32	0.5
7KG-G019-11	ZB-XLB	60	0.25	0.25
7EG-G024-02	ZB-XLB-HT	15	0.25	0.1
7EG-G024-11	ZB-XLB-HT	15	0.25	0.25
7FD-G024-08	ZB-XLB-HT	20	0.18	0.18
7HG-G024-02	ZB-XLB-HT	30	0.25	0.1
7HG-G024-11	ZB-XLB-HT	30	0.25	0.25
7HM-G024-02	ZB-XLB-HT	30	0.32	0.1
7HM-G024-11	ZB-XLB-HT	30	0.32	0.25
7KG-G024-11	ZB-XLB-HT	60	0.25	0.25
7KG-G045-10	Zebron ZB-Dioxin	60	0.25	0.2
7PD-G045-47	Zebron ZB-Dioxin	40	0.18	0.14
7FD-G044-47	Zebron ZB-PAH-CT	20	0.18	0.14
7HG-G044-10	Zebron ZB-PAH-CT	30	0.25	0.2
7PD-G044-47	Zebron ZB-PAH-CT	40	0.18	0.14
7AB-G000-00-GZ0	Zebron Z-Guard	5	0.1	0
7AD-G000-00-GH0	Zebron Z-Guard	5	0.18	0
7AD-G000-00-GZ0	Zebron Z-Guard	5	0.18	0
7AD-G000-00-GZK	Zebron Z-Guard	5	0.18	0
7AG-G000-00-GH0	Zebron Z-Guard	5	0.25	0
7AG-G000-00-GHK	Zebron Z-Guard	5	0.25	0
7AG-G000-00-GZ0	Zebron Z-Guard	5	0.25	0
7AG-G000-00-GZK	Zebron Z-Guard	5	0.25	0
7AK-G000-00-GH0	Zebron Z-Guard	5	0.53	0
7AK-G000-00-GHK	Zebron Z-Guard	5	0.53	0
7AK-G000-00-GM0	Zebron Z-Guard	5	0.53	0
7AK-G000-00-GZ0	Zebron Z-Guard	5	0.53	0
7AK-G000-00-GZ1	Zebron Z-Guard	5	0.53	0
7AK-G000-00-GZK	Zebron Z-Guard	5	0.53	0
7AM-G000-00-GH0	Zebron Z-Guard	5	0.32	0
7AM-G000-00-GHK	Zebron Z-Guard	5	0.32	0
7AM-G000-00-GZ0	Zebron Z-Guard	5	0.32	0
7AM-G000-00-GZK	Zebron Z-Guard	5	0.32	0
7CB-G000-00-GZ0	Zebron Z-Guard	10	0.1	0
7CD-G000-00-GZ0	Zebron Z-Guard	10	0.18	0
7CG-G000-00-GH0	Zebron Z-Guard	10	0.25	0
7CG-G000-00-GHK	Zebron Z-Guard	10	0.25	0
7CG-G000-00-GZ0	Zebron Z-Guard	10	0.25	0
7CG-G000-00-GZK	Zebron Z-Guard	10	0.25	0
7CK-G000-00-GH0	Zebron Z-Guard	10	0.53	0
7CK-G000-00-GHK	Zebron Z-Guard	10	0.53	0
7CK-G000-00-GZ0	Zebron Z-Guard	10	0.53	0
7CK-G000-00-GZK	Zebron Z-Guard	10	0.53	0

Наименование	Фаза	Длина (м)	Внутренний диаметр (мм)	Толщина пленки (мкм)
7CM-G000-00-GH0	Zebron Z-Guard	10	0.32	0
7CM-G000-00-GZ0	Zebron Z-Guard	10	0.32	0
7CM-G000-00-GZK	Zebron Z-Guard	10	0.32	0
7JG-G000-00-GZ0	Zebron Z-Guard	50	0.25	0
7JK-G000-00-GZ0	Zebron Z-Guard	50	0.53	0
7JM-G000-00-GZ0	Zebron Z-Guard	50	0.32	0
KG0-7868	Zebron Z-Guard	2	0.25	0

Архангельск (8182)63-90-72
Астана (7172)727-132
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13

Сургут (3462)77-98-35
Тверь (4822)63-31-35
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Ярославль (4852)69-52-93

www.phenomenex.nt-rt.ru || pxp@nt-rt.ru